

# PICASSO LITE+®

Picasso Lite+® is the most celebrated, endorsed, and awarded dental laser in the industry. Its ease of use and affordable price point make it the perfect solution for treating soft tissue. The Picasso Lite+ laser was designed with the first-time laser dentist or hygienist in mind while simultaneously offering the necessary capabilities for experienced practitioners looking to create masterpiece smiles and unleash their inner artist. Picasso Lite+ laser is an effective and affordable solution for every operator in a dental practice.

## FEATURES

- Easy to Use Presets
- Wireless Foot Control
- Built-In Treatment Timers for Periodontal Treatment and Aphthous Ulcer Treatment
- Multiple Languages Menu
- Continuous/ Pulsed Modes

## APPLICABLE PROCEDURES

- Frenectomy
- Fibroma Removal
- Troughing
- Hemostasis
- Periimplantitis
- Tissue Tag
- Coagulation
- Sacular Debridement
- Curettage
- Lip and Tongue Tie
- Herpetic Lesions
- And More...



VISIBLE AND INVISIBLE LASER RADIATION. AVOID EYE OR SKIN EXPOSURE TO DIRECT OR SCATTERED RADIATION.  
 WAVELENGTH: 810NM WAVELENGTH: 650NM  
 OUTPUT: 50W OUTPUT: 50W  
 CLASS 4 LASER PRODUCT CLASS 2 LASER PRODUCT  
 RAYONNEMENT LASER VISIBLE ET INVISIBLE. EXPOSITION DANGEREUSE DE L'ŒIL OU DE LA PEAU AU RAYONNEMENT DIRECT OU DIFFUS  
 LA LONGUEUR D'ONDE: 810NM APPAREIL A LASER DE CLASSE 4  
 ÉMISSION MAXIMALE: 50W APPAREIL A LASER DE CLASSE 2  
 LA LONGUEUR D'ONDE: 650NM APPAREIL A LASER DE CLASSE 2  
 ÉMISSION MAXIMALE: 50W  
 IEC 60825-1:2014 - SAFETY OF LASER PRODUCTS  
 REV. 1

## TECHNICAL SPECIFICATIONS

WEIGHT	DIMENSIONS	OUTPUT VOLTAGE	OUTPUT CURRENT	OUTPUT POWER RANGE	LASER CERTIFICATION	MEDIUM	WAVELENGTH	POWER MODES	DELIVERY FIBER	OUTPUT RANGE
2.0 lbs (1 kg)	9.5" x 6" x 6.2"	+12.0 V DC	5.0 AMPS	60 Watts MAX	IV (4)	GaAlAs	810 ± 10nm	Continuous, Pulse Repeat	200µm - 400µm	0.5-3.0 Watts ±20% output



# THE MOST AWARDED DENTAL LASER IN THE INDUSTRY

The Picasso Lite+ Diode Laser was designed for versatility, usability, precision, and affordability. Offering a wide range of necessary capabilities for soft tissue surgical procedures and periodontal treatments, the Picasso Lite+ soft tissue laser is ideal for both dentists and hygienists alike. The Picasso Lite+ is equipped with a multi-tip handpiece, wireless foot control, and customizable modes & presets—and all at the most affordable price the market today. ***We guarantee it!***



# BETTER HEALING IS BETTER RETENTION

Diode lasers help reduce tissue damage, bleeding, and receding gums. The Picasso Lite+ Laser streamlines the treatment of periodontal disease and minimizes the need for local anesthetic, stitches, and periodontal packs.

Using soft-tissue diode lasers for surgical procedures in your dental practice increases revenue and retention by decreasing time in the dentist chair and expediting patient healing.

The Picasso Lite+ Diode Laser can perform up to 32 different periodontal and soft-tissue laser procedures, allowing you to expand the services you offer. From gingivectomies to frenectomies, never outsource periodontal laser procedures again.

