

SHEER/WHITE!®

PROFESSIONAL
TEETH WHITENING SYSTEM
Available ONLY From Your Dentist

NMD
INNOVATING
DENTISTRY
SINCE 2000



whiten anytime anywhere



Apply



Relax



Drink



Exercise



Talk

Sheer White!® professional whitening strips stick and mold perfectly to the teeth, becoming almost invisible. They're comfortable to wear during any activity and cause little-to-no-sensitivity.

SHEER WHITE!®

**PROFESSIONAL
TEETH WHITENING SYSTEM**
Available **ONLY** From Your Dentist

Offer your patients a powerful teeth whitening solution that is fast, effective, and virtually painless. *Sheer White!* allows you to offer whitening to all your patients by giving you an affordable option compared to other professional whitening products. Custom trays and gels can be painful, and too pricey for some of your patients. *Sheer White!* will save chair time, save money, and help create repeat whitening visits. Your patients will love their whiter smile as much as they love saving time and money. ***Better Retention, Better Revenue, Better Whitening!***

In-Office



Apply an In-Office whitening strip at the end of any appointment; application only takes one minute. The *Sheer White!* In-Office strip is a 20% Hydrogen Peroxide formula delivered by our patented non-slip strip, and works with minimal sensitivity. *Sheer White!* In-Office is capable of whitening up to 3 shades whiter in just 30-60 minutes.

6 VALUE INTRO KIT 006-00251

24 PATIENT BULK PACK 006-00293

Take-Home



Send your patients home with the *Sheer White!* Take-Home patient treatment kit, a 20% carbamide peroxide formula delivered by our patented go-anywhere non-slip whitening strip. Each *Sheer White!* Take-Home kit includes a 5-day treatment, consisting of 10 strips. Patients' teeth will whiten up to 5 shades whiter in just 5 days.

6 PATIENT INTRO KIT 006-00081

24 PATIENT BULK PACK 006-00294

For best results, combine the *Sheer White!* In-Office and Take-Home strips and use as one simple complete whitening system:

- Apply the In-Office strips at the end of any appointment. Patients wear for 30-60 minutes for whitening up to 3 shades.
- Give your patients a *Sheer White!* take-home kit and have them follow instructions for up to 5 more shades in 5 days.